

# **PsyClick: Clinical Decision Support System for Detecting Single-Session Psychomotor Disturbances via EWMA-Based Multivariate Analysis of Keystroke and Cursor Dynamics**

by

Ayaay, Martonie Loreign L.  
Ballano, Denise Nicole M.  
Huertas, Lana Denise R.  
Tablate, Judea Ann C.

## **Technical Documentation**

FEU Institute of Technology  
July 2026

©2026 Ayaay, Martonie Loreign L., Ballano, Denise Nicole M.,  
Huertas, Lana Denise R., Tablate, Judea Ann C.  
All Rights Reserved

### **1. Purpose and Scope**


This document is the consolidated Technical Documentation for PsyClick Secure, required under the CS0029 Final Project Specification (Section IV.B, item 3). It contains the System Architecture Diagram, Database Schema, API Documentation, Installation Guide, and a condensed User Manual for the psychlick-secure application.

The full clinician and patient operating guide is provided as a separate User Manual deliverable; the section below summarizes the operating procedure.

### **2. System Architecture Diagram**

PsyClick Secure applies defense in depth within a single-workstation trust boundary. A local Electron/React presentation layer communicates with a loopback Flask API (127.0.0.1:5101); server middleware authenticates bearer sessions; role checks and clinician ownership scope protected records; a fail-closed consent gate precedes telemetry capture; SQLite persistence uses parameterized queries and transactions; audit events are hash-chained; and Windows DPAPI protects exports and backups.

**Figure 38. PsyClick Secure — Layered Defense-in-Depth Architecture**



*Figure 1. PsyClick Secure layered security architecture.*

**3. Database Schema**


The secure database classifies clinician identity, sessions, consent, clinical/behavioral records, audit evidence, and recovery metadata. Direct SQL string concatenation is prohibited; parameterized queries, foreign-key constraints, transaction boundaries, and application-level authorization are required.

| Entity            | Primary Key  | Security-Relevant Fields                                                                | Classification                |
|-------------------|--------------|-----------------------------------------------------------------------------------------|-------------------------------|
| clinicians        | clinician_id | script password_hash;<br>role; status;<br>failed_attempts;<br>locked_until              | Restricted identity           |
| security_sessions | session_id   | token_hash; clinician_id;<br>role; last_seen_at;<br>expires_at; revoked_at              | Confidential session metadata |
| consent_records   | consent_id   | patient_id; clinician_id;<br>consent_version;<br>decision; recorded_at;<br>withdrawn_at | Restricted consent            |

|                 |            |                                                                          |                                   |
|-----------------|------------|--------------------------------------------------------------------------|-----------------------------------|
| intake_sessions | session_id | patient/session identifiers;<br>features; scores; flag;<br>rationale     | Restricted<br>clinical/behavioral |
| audit_log       | log_id     | actor; actor_id; action;<br>outcome; prev_hash;<br>entry_hash; timestamp | Confidential<br>accountability    |
| backup_records  | backup_id  | created_by; file_path;<br>SHA-256; encrypted;<br>status; verified_at     | Confidential recovery<br>metadata |

The security-aware entity-relationship diagram is shown below.

**Figure 40. PsyClick Secure — Core Entity-Relationship Diagram**



*Figure 2. Core entity-relationship diagram.*

## 4. API Documentation

All endpoints are served from the loopback Flask API at <http://127.0.0.1:5101/api/> and require an authenticated bearer session except for login, first-account bootstrap registration, and the readiness ping. Access is enforced by role-based middleware.

### Authentication

POST /api/login, POST /api/register (bootstrap or admin-invoked), POST /api/verify-clinician, POST /api/logout

### Clinical workflow

POST /api/intake/start, GET /api/next-client-id, GET /api/next-tester-id, calibration (keyboard/mouse start and save), PHQ-9/GAD-7 and emotional-task assessment routes, POST /api/assessment/finish

### Records and reporting

GET /api/stats, GET /api/sessions/recent, GET /api/session/<id>, GET /api/patients, GET /api/patients/<id>/sessions, POST /api/export/report, POST /api/export/summary

### Audit and security monitoring

GET /api/audit (admin, auditor), POST /api/audit/log, GET /api/audit/verify (admin, auditor) — hash-chain verification, GET /api/security/intrusion-alerts (admin, auditor)

### Administration

GET /api/admin/users (admin), PATCH /api/admin/users/<id> (admin), POST /api/admin/backup, POST /api/admin/backup/verify, GET /api/admin/security-status (admin)

### Sync and health

GET /api/ping, GET /api/db-health, GET /api/sync/status, POST /api/sync/now, POST /api/sync/pull (optional Supabase cloud sync; no-op if unconfigured)

### Normative / research mode

POST /api/normative/login, POST /api/normative/mode, GET /api/normative/stats, POST /api/normative/compute, GET /api/normative/compare/<id>

## 5. Installation Guide

On a supported Windows workstation:

- Create a Python virtual environment and install requirements.txt.
- For a packaged/controlled deployment, set the PSYCLICK\_SECURE\_HOME environment variable for the secure data directory (defaults to %APPDATA%\PsyClickSecure).
- Build the backend and frontend with BUILD\_SECURE.bat (PyInstaller creates dist-python\psychlick\_api\psychlick\_api.exe; Vite builds the frontend).
- To run in development, start api\_server.py on 127.0.0.1:5101, install frontend dependencies, and run npm run dev — or use start\_secure.bat.
- To produce the installable desktop app, run npm run dist:installer in the frontend folder.
- Launch the application and create the first administrator account; further accounts require an administrator session.

*The separate source archive and frontend build archive are supplied. A fresh independent-workstation installation still requires recorded evidence.*

## 6. User Manual (Summary)

The full Clinician Manual and Patient Guide are provided as a separate User Manual deliverable. The core operating procedure is:

- Sign in with an individual clinician account and confirm the API readiness indicator.
- Create a pseudonymous client and record informed consent (capture is blocked until consent is recorded).
- Run keyboard and mouse calibration, then administer PHQ-9, GAD-7, and the emotional-task prompts.
- Review the dashboard — T<sup>2</sup>, PSI/PAI, confidence, and the Green/Amber/Red signal — alongside clinical observation.
- Export only to an approved protected destination; exports and backups are DPAPI-encrypted with a SHA-256 check.
- Review the audit log, log out, and lock the workstation.

## 7. Approval

| Name | Signature |
|------|-----------|
|------|-----------|

|                            |  |
|----------------------------|--|
| Ayaay, Martonie Loreign L. |  |
| Ballano, Denise Nicole M.  |  |
| Huertas, Lana Denise R.    |  |
| Tablate, Judea Ann C.      |  |